

NEW

Bench Top Supreme Recycler

*Twice the Recycler.
Half the Space.*

*Designed for use
on the counter!*

*Fast and more
automatic!*

*Does the work of a
floor model!*

**Recycle and Reuse
the following:**

- *Xylene*
- *Xylene Substitutes*
- *Alcohol*
- *Formalin*



CBG BIOTECH

Capacity

- 2.5 gallons

Footprint

- 12.75"w x 23"d x 39"h

Cycle Times

- Xylene – 1 hour per gallon
- Alcohol – 1.5 hours per gallon
- Formalin – 2-2.5 hours per gallon

Unit Features

- Fits conveniently on counter
- All front loaded
- Built-in spill containment
- Automatic waste drain
- Air cooled – no ventilation needed
- No glass parts
- No water hookup
- No hot surfaces
- Listed to UL 2208

Unit Specifications

- Fractional distillation
- Electrical hook-ups: 120vac, 50/60 HZ, 10-15 amps
- Satisfies the latest uniform fire code requirements

Responsive Customer Service

- 30 day free trial, with no hidden costs
- Full recycler replacement and parts while under warranty or service contract
- Technical staff available 24 hours a-day, 7 days a-week
- Toll free service number available at 1-800-941-9484

Pure &
Simple.



FREE TRIAL!

For a no obligation 30 day evaluation, call 1-800-941-9484. See for yourself how CBG Biotech can revolutionize your laboratory.

CBG BIOTECH

2211 Lake Club Drive, Suite 205
Columbus, Ohio USA 43232
614-863-1650 • 1-800-941-9484
fax 614-863-1676
e-mail info@cbgbiotech.com
www.cbgbiotech.com

Tested to UL2208 safety and electrical specifications for solvent distillation units. The Bench Top Supreme is the smallest unit today that's been tested to these regulations.