



Recycle and Reuse – Green, Safe and Easy!

New! 4-LITER COMPACT SOLVENT RECYLER

- RECYCLE and RE-USE Xylene, Alcohol and other solvents used in the Lab!
- Weighs less than 30 kgs.
- ONLY 36 cm wide x 40.5 cm deep x 82 cm high!
- Fits conveniently on counter
- Minimal Operator Time – Fully Automated!
- Air cooled – no ventilation needed
- No glass parts or water hookup
- Built to CE Standards!
- No Hot Surfaces!



Approx. Cycle Times (full tank): Xylene – 2 hours Alcohol – 3 hours

High Purity and Recovery Rates!

Even small labs can save money,

Get great results and

GO GREEN!