



Free Trial!

*Ask for a no-obligation evaluation
at your laboratory. Revolutionize your
Laboratory to Get Pure Profit!*

CBG BIOTECH

2211 Lake Club Drive, Suite 205
Columbus, Ohio USA 43232
800-941-9484 • 614-863-1650
fax 614-863-1676
e-mail info@cbgbiotech.com
www.cbgbiotech.com

CBG BIOTECH

Solvent and Formalin Recycling



Pure & Simple.



Solvent and Formalin Recycling

Recycle and Reuse at least 90 to 95% of the following:

- Xylene
- Xylene Substitutes
- Alcohol
- Acetone
- Formalin

CBG Biotech Recyclers are designed and built to offer the **highest purity**, while also being **safe** and **easy to use**. By reducing your solvent purchases and your waste disposal, handling and storage costs, you can start realizing significant savings IMMEDIATELY! We provide our customers with economical ways to achieve a **cleaner environment** and **safe working areas**.

Pure & Simple

5 Gallon FormaSolve



10 Gallon Solvent



2 1/2 Gallon Solvent



Compact



CBG BIOTECH

Listed to UL Standards 2208 & 3101, CSA (Canadian Standards), and CE (European Marking) by ETL

Unit Specifications

- Dimensions:
 - 2 1/2 gallon Recycler – 24" wide, 67" high, 30.5" deep
 - 5 gallon Recycler – 28.5" wide, 67" high, 30.5" deep
 - 10 gallon Recycler – 31" wide, 67" high, 30.5" deep
 - FormaSolve Recycler (2 1/2 or 5 gallon) – 28.5" wide, 67" high, 30.5" deep
- Type of distillation - fractional
- Type of control: microprocessor, fully automated
- Electrical hook-ups: 110-120v., 60 cycles, 7-15 amps
- Satisfies the latest Uniform Fire Code requirements for solvent distillation units
- Listed to UL Standards 2208 & 3101, CSA (Canadian Standards), and CE (European Marking) by ETL

Unit Features

- Automatic waste disposal (except compact unit)
- Automatic decontamination process at solvent change
- Automatic start
- No cool-down time required between runs
- No water hook-up or special ventilation required
- No need to change lines or pots
- No hot surfaces – Cool to the touch
- No glass parts
- Mobile – unit has locking wheels (except compact unit)
- Electronic fill pump on 5 gallon or larger capacity units
- Built-in spill containment
- Variable volume selection

Solvent/Formalin	Cycle Time	Purity	Recovery
Xylene/Xylene Substitutes	1 hour per gallon	>99%	95%
Alcohols	1.5 hours per gallon	99% <small>(Concentration 92-99% depending on collection)</small>	90%
Formalin	2 hours per gallon	99% <small>(Concentration 10-12%)</small>	80-90%

Responsive Customer Service

- Board certified histology and laboratory professionals to provide knowledgeable support
- 24 hour 1-800 customer service telephone support
- On-site service available
- Warranty replacement packages available

Options

- Free 30 Day Trial
- Compact Benchtop Models available
- Customization available to meet specific needs
- Vapor Detector
- Fire Extinguisher
- Drum Dollies for easy transport
- Solvent Self-Selection Program

Financing Options

Start saving immediately! CBG Biotech will work with you to provide short-term financing options to overcome short-term capital expenditure and budget constraints. CBG Biotech also offers long-term lease financing, including lease-to-own options. Leases are available for as short as one year, and for as long as 5 or 6 years, with or without an option to purchase the equipment at the end of the lease period. Full service support is provided as part of the lease arrangement.

Now Distributing CBG Biotech Recyclers Worldwide!