

NEW

Thermo Scientific PrintMate ES Cassette Printing System


Meeting the needs of lean workflow and error reduction within your laboratory

Eliminate waste and streamline your labeling process by printing cassettes on-demand with the PrintMate ES cassette printer at the grossing station.

- Prints standard cassettes and **NEW** Thermo Scientific SecureSette tissue cassettes
- Small footprint allows the printer to fit easily next to the grossing station
- Eliminate manual handwriting processes
- Print only the cassettes you need, right at the point of use
- Print barcodes on cassettes to automate downstream processes



NEW INNOVATIVE TISSUE CASSETTE



Enable lean workflow and enhance patient safety with the NEW SecureSette™ tissue cassette

The SecureSette tissue cassette allows final tissue orientation during the gross dissection process

- Tissue sample remains oriented inside the cassette through embedding which improves sample chain of custody.
- The need for tissue orientation at embedding is eliminated enhancing workflow.
- Embedding time is reduced by 50%.

Thermo Scientific PrintMate ES Cassette Printer

Order Information

Description	Order Number
PrintMate ES 150 Cassette Printer	A82800001
PrintMate ES 450 Cassette Printer	A82800002
PrintMate ES 900 Cassette Printer	A82800003

Includes PrintMate ES, USB cable, power cable, 2, 6 or 12 standard cassette drawers, one SecureSette drawer (two for PrintMate ES 900), one roll of tape (installed), LabWriter software CD, operator guide

Description	Order Number
PrintMate ES 150 Cassette Printing System	A82800004
PrintMate ES 450 Cassette Printing System	A82800005
PrintMate ES 900 Cassette Printing System	A82800006

Includes PrintMate ES, USB cable, power cable, 2, 6 or 12 standard cassette drawers, one SecureSette drawer (two for PrintMate ES 900), one roll of tape (installed), LabWriter software CD, operator guide, All-in-one PC, digital imager

PrintMate ES Accessories

Description	Order Number
Collection System (includes cable, 4 trays for standard cassettes, 4 trays for SecureSette and output guides for both types of cassettes)	A79510014
Collection trays, standard cassettes (box of 4)	A79510022
Collection trays, SecureSette (box of 4)*	A82810020
Hot Foil Tape, box of 5 rolls	A79510013
All-in-one Touch screen PC	A79510021
Digital Imager (Linear/2D Barcode Reader)	B35310006
Cassette hopper, standard cassettes (ea)	A79510024
Cassette hopper, SecureSette (ea)	A82810021
Collection system output guide, SecureSette* (2 ea)	A82810022
Collection system output guide, standard cassettes (2 ea)	A82810023
Cassette guide chute with extension, SecureSette* (2 ea)	A82810024

* SecureSette specific accessories also accommodate standard cassettes



Specifications

Dimensions: (DxWxH):

PrintMate ES 900 w/o cassette tubes	15 x 17 x 11.4 in. (38 x 43.5 x 29cm)
PrintMate ES 900 w/ cassette tubes	15 x 17 x 32 in. (38 x 43.5 x 81cm)

Weight:

PrintMate ES 900 w/o cassette tubes	31 lb (14 kg)
PrintMate ES 900 w/ cassette tubes	37.5 lb (17 kg)

Power Requirements:

100 – 240 VAC, 50/60 Hz,
150 VA (512 BTU/h)

Print Media Type:

Thermal transfer ribbon

Print Media Capacity:

Approximately 20,000
cassettes per roll

Print Media Color:

Black

Throughput:

12 seconds per
cassette typical

Cassette Capacity:

PrintMate ES 150	150 cassettes
PrintMate ES 450	450 cassettes
PrintMate ES 900	900 cassettes

To learn more about how on-demand cassette printing can impact in your laboratory, call: (800) 522-7270 or visit www.thermoscientific.com/printmate