

# Thermo Scientific SlideMate Slide Printer

## Bringing On-Demand Slide Printing to the Microtome!

Studies show that thousands of specimen identification errors occur annually which may lead to adverse patient events. A large percentage of these errors are due to improperly labelled slides. The small SlideMate™ footprint (about the size of a toaster) brings the convenience of on-demand slide printing directly to the microtome. This drastically eliminates errors by allowing the laboratory to reliably print accurate information at the moment the slide is required.

- On-demand slide printing reduces errors
- Prints directly onto the slide eliminating the need for labels or handwriting
- Small footprint fits easily next to the microtome
- Chemical resistant media ensures print reliability
- Flexible data inputs allow easy LIS integration



***“Slide labeling alone accounted for 67% of defects”***

Source: The Henry Ford Production Systems Defects and Waste in Surgical Pathology as a Basis for Quality Improvement Initiatives

©2007 American Society for Clinical Pathology

## SlideMate™ Offerings for your Laboratory

- On-demand slide printing streamlines workflow and increases accuracy of identification
- Small footprint fits easily into the lab
- Flexible data input methods make for easy integration into the laboratory
- High-quality printing resistant to harsh chemicals
- Slides are held inside the unit's side door to protect them from dust and the environment
- User can select the number of slides to print per case, which enhances productivity in the lab

### Specifications

Width	6.0", 152mm
Height	8.9", 226mm,
Depth	10.4", 264mm,
Weight	4.3kg, 9.5lbs.
Power Requirements	Universal 100-240 VAC 47-63 Hz 1.2A Maximum
Color	Sterile Steel
Slide Type	3"x1" or 25mm x 75mm x 1mm with Superfrost or Colorfrost end
Slide Capacity	150 pieces (input stack) 10 pieces (output tray)
Media Type	Thermal transfer ribbon, 1" wide roll
Media Capacity	1 roll, approximately 12,000 Slides
Media Color	Black
Throughput	3-6 Slides/minute typical

Streamline workflow and minimize specimen identification errors in your laboratory with the Thermo Scientific PrintMate cassette printer.



### Order Information

Description	Order Number
<b>SlideMate Print On-Demand Printing System</b> includes SlideMate Printer, Digital Imager, Keyboard Interface, Mini Keyboard, Box of Slides	<b>B81300001</b>
<b>SlideMate with LabWriter™ Software</b> includes SlideMate Printer, LabWriter Basic Software for SlideMate, Box of Slides	<b>B81300002</b>
<b>SlideMate with PC – Premium Package</b> includes SlideMate Printer, All-in-one PC, Digital Imager, LabWriter Basic Software, Box of Slides	<b>B81300003</b>
<b>SlideMate Printer</b> includes SlideMate Printer	<b>B81300004</b>
SlideMate Available Accessories	Order Number
Digital Imager	B35310006
All-in-one PC	A79510021
LabWriter Basic Software for SlideMate	B81310001
Network Cable	B35310012
Crossover Ethernet Cable	B81310005
Keyboard Interface	B81320043
Mini Keyboard	B81310004

© 2009 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

**Anatomical Pathology** 4481 Campus Drive  
Kalamazoo, MI 49008  
(800) 522-7270

Robert-Bosch-Str. 49  
D-69190 Walldorf/Germany  
+ 49(0)6227-8360

93/96 Chadwick Road  
Astmoor, Runcorn Cheshire WA7 1 PR  
0800 018 9396, + 1 (0) 1928 562501

[www.thermo.com/pathology](http://www.thermo.com/pathology)